

1 ACLS SQL8.5.x and SQL8.6.x bug-fixes and new features:  
2  
3 1. Resource MoxyTap manager  
4 - New module to work with moxyTap NFC card tap device for sign in and sign out to  
the resource and lab space  
5 2. Lab title styling customised color in global setting:  
6 - LabTitleBackgroundColor  
7 - LabTitleFontColor  
8 3. Message/Alert/Pinboard icon  
9 - Restyling in blue balloons  
10 4. Mobile app  
11 - Shown resource website link as active link other than text on resource page  
12 - Updated user access date when accessing app  
13 5. User profile manager  
14 - Added registered date to the search table  
15 - Renamed "created" as "registered"  
16 - Fixed the glitch for the account selection to the account name using speech quotes  
17 6. Sample tracking  
18 - Fixed the glitch to set up sample operation groups for sample-based and time-based  
selection  
19 7. Resource trainer on Dashboard  
20 - Added resource location to the trainer table  
21 8. Service performance enhancement  
22 - Improve the scripts for system daily cache management  
23 9. UWA RIMS  
24 - Modification per RIMS 2024 proposal  
25 10. PHP engine  
26 - Updated to the latest v8.4.1  
27 11. Enterprise service - Business intelligence (BI)  
28 - Added feature: crossed-ref reports cover the booking usages and onsite tracker  
usages, and including the comparsion between bookings and onsite tracker  
sign-in/sign-out.  
29 - Improved feature: issued user certificates including the history records for every  
users  
30  
31  
32  
33  
34  
35  
36