

1 ACLS SQL8.7.x bug-fixes and new features:  
2  
3 1. Cybersecurity features  
4 - Added OTC pages to online reg, password reset, access denied pages  
5 - Added new anti-hacking and anti-scanning for login page of non-enterprise ACLS  
6 instance  
7 2. moxyTap features  
8 - Working with NFC card reader device to track the real-time usages of the selected  
9 resources  
10 - NFC card reader devices provided by UniLab, tap to sign in and tap to sign out,  
11 supporting lab spave sign in and single instrument sign in  
12 3. Terms and Conditions file support  
13 - Discontinue text file support  
14 - Added PDF viewer for terms and conditions  
15 4. Access group manager  
16 - Added deafult user group setting for user reg approval  
17 5. UWA RIMS  
18 - Fixed the bugs as per modification per RIMS 2023 proposal  
19 6. MCN  
20 - Fixed the bugs as per modification per MCN 2020 proposal  
21  
22  
23  
24